
9/96-Gan/Sf-

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This apparatus, consisting of the light velocity measuring instrument and 2 triple mirrors, permits measurement of the velocity of light from the path and the transit time of extremely short light pulses. These are reflected in a large triple mirror so that they pass through the measurement arrangement twice and are then converted to voltage pulses, which are displayed on the oscilloscope for determination of the transit time.

This configuration uses measuring distances of 10 m to 20 m. To permit measurement of the transit times of approx. 65 ns to 130 ns using standard school oscilloscopes, the light velocity measuring instrument is equipped with a trigger output and a quartz-stabilized 10-MHz output. The signals for external triggering of the oscilloscope are optimally matched to the voltage pulses, so that complete pulses can be displayed at all times. Where necessary, the signal of the 10-MHz output can be used on the second oscilloscope channel as a calibrated time base; the edges of the signal (period 100 ns) then serve as a measuring standard which is independent of the properties of the individual oscilloscope.

Thanks to their narrow peak width of 20 ns, the voltage pulses obtained from the light pulses can also be used for experiments on the propagation of signals in coaxial cables, e.g. determining the propagation velocity ( $2/3$  of the velocity of light); reflection at an open and a terminated end; the effect of a terminating resistance; determining the wave resistance.

#### Literature:

Experiment descriptions in "Leybold Physics Leaflets", Part 2 (597 522)

### 1 Safety notes

- Protect the light velocity measuring instrument from mechanical shocks (e.g. knocks or falls) to prevent disturbing the optical adjustment of the transmitter LED, the beam divider and the receiver diode.
- Protect all apparatus from dust and smoke to prevent reduced light yield due to dirty optical elements.
- Store the light velocity measuring instrument upright (e.g. in a stand base, 300 02); protect opening (1.7) with cover (1.8).
- Protect the triple mirror against scratching or similar damage.
- Connect only  $50\text{-}\Omega$  coaxial cable to the outputs (output resistance  $50 \Omega$ ) to ensure distortion-free display of the signal.

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## EMC information

The light velocity measuring instrument is a class A device within the meaning of the standard EN 55011. It can cause radio interference in residential areas. In this case, the operator may be obligated to implement appropriate remedial measures at his or her own expense.

When operating the device in an experiment room in general or vocational schools or other educational institutions, however, it may be assumed that in general no interference will occur in residential areas, providing that a protective distance of 30 m is maintained and the device is operated for only brief periods.

*Measures for ensuring compliance with the requirements of class A according to standard EN 55011:*

- Use only screened measuring leads of the shortest possible length, e.g. BNC cable, 1 m (501 02).
- Do not switch on the light velocity measuring instrument ① by connecting the plug-in power supply until you have carried out all experiment preparations which do not require a light beam.
- Do not operate the device longer than required for carrying out the experiment.

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## 2 Measuring principle

### 2.1 velocity of light (see Fig. 1)

A high-performance LED (a) driven by an oscillator circuit emits short red light pulses with a repetition frequency of 40 kHz. The beam passes through a semitransparent mirror which acts as a beam divider (b); one half is reflected as a "reference beam" (see below) toward housing window (c<sub>0</sub>); the other half, the "measuring beam", passes through window (c<sub>1</sub>). A lens (f = 200 mm) images the LED on a triple mirror (d<sub>1</sub>), which is (by definition) at a distance s/2; this mirror reflects the beam precisely into itself and, via beam divider (b), into a receiver diode (e). This converts each light pulse that has travelled the path s into a voltage pulse, which is displayed on the oscilloscope as the "time signal" U<sub>1</sub>.

If the triple mirror is moved by  $\Delta s/2$ , so that the light path changes by  $\Delta s$ , the time signal U<sub>1</sub> on the oscilloscope is shifted by  $\Delta t$ . The velocity of light can be determined from the slope of the line  $\Delta s = f(\Delta t)$  by recording multiple value pairs.


To measure the absolute time t which a light pulse requires to travel path s, a reference signal U<sub>0</sub> for the zero point is required in addition to the time signal U<sub>1</sub>; this reference signal is generated using a second, small triple mirror (d<sub>2</sub>). This is set up either on window (c<sub>0</sub>) in the light beam which is reflected upward by the beam divider (b) or directly in the measuring beam at win-

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dow ( $c_1$ ), and reflects the reference beam to receiver diode ⑥. The two light paths are equivalent. Placing the mirror in the measuring beam is recommended for didactic reasons; the inevitable impairment of the measuring beam is obviated when window ( $c_0$ ) is used. The reference beam which is converted into a voltage by the receiver diode is displayed on the oscilloscope as reference signal  $U_0$ . The time difference between  $U_1$  and  $U_0$  is the time  $t$  which the light requires to travel path  $s$ .

## 2.2 Transit time of signals in coaxial cables

When studying transit times in coaxial cables, the signals generated in window ( $c_0$ ) by reflection of the light pulses at the small triple mirror ( $d_2$ ) are fed to the oscilloscope via a T-connector as the reference signals  $U_0$ , and are also fed into a coaxial cable with the length  $l$  ( $l \geq 10$  m) as the measuring signals  $U_1$ . The  $U_1$  signals pass through the cable and are reflected at the open end without phase shift or – when the inner and outer conductors are short-circuited or connected via a resistor – reflected at the terminated end with phase shift. After travelling path  $2l$ , these reach the oscilloscope shifted by transit time  $t$  with respect to the  $U_0$  signals.

## 3 Scope of supply, description, technical data (Fig. 2)

Fig. 2

Lichtgeschwindigkeits-Meßgerät (476 50)  
Light velocity measuring instrument (476 50)

### ① Light velocity measuring instrument

This apparatus comprises

- the pre-adjusted optical elements transmitter LED, beam divider and receiver diode (see also section 2) plus
- the electronics for generating the pulses available at outputs (1.1) to (1.3) (output resistance  $50 \Omega$ ) (see Fig. 3/4).

#### (1.1) Pulse output

Supplies the voltage pulses obtained from the transmitter LED via the receiver diode as time signals for time measurement on the oscilloscope (see Fig. 3).

Pulses:

- Repetition frequency: 40 kHz
  - Peak width: 20 ns
  - Half-value peak width: typically 5 ns
- Wavelength of light: 615 nm

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Phase relation of signals at the pulse and trigger outputs (schematic representation)

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Spiegelfläche; 20 mm x 35 mm (1 Element)

### (1.2) Trigger output

The pulses trigger the oscilloscope externally. The positive edge of the square-wave trigger voltage (TTL-signal) leads the emission of the light pulses by about 60 ns (see Fig. 3). As a result, the complete measuring signals can be displayed on standard school-quality oscilloscopes even for extremely short light paths (e.g. triple mirror directly in front of windows (1.6) or (1.7)).

### (1.3) 10-MHz output

The quartz-stabilized AC voltage signals with a period of 100 ns can be used on the second oscilloscope channel as a calibrated time base. The edges of the signal serve as a measurement standard:

- when the time measurement in the calibrated ranges of the oscilloscope cannot be carried out under optimum use of the oscilloscope screen or
- when the time base of the oscilloscope is to be checked and/or corrected.

The signal shape depends on the bandwidth of the oscilloscope used. For a 30 MHz oscilloscope, this approximately corresponds to Fig. 4; for a 100 Hz instrument the edges are steeper, the signals resemble a square-wave.


Fig. 4

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### (1.4) Phase adjuster

Permits a phase-shift of the 10-MHz signal over 1 period for matching to the phase relation of the measuring signal (optimizing the evaluation conditions, see section 4.3.3).

### (1.5) Switched socket for connecting the plug-in power supply with cannon plug 12 V, 20 W (562 791 for 230 V AC or 562 792 for 115 V DC)

### (1.6) Window for measuring beam

### (1.7) Window for reference beam

### (1.8) Cover for window (1.7)

Housing dimensions: 115 mm x 65 mm x 155 mm

Stand rod: 30 mm long; 10 mm dia.

Weight: 1.4 kg

## ② Pair of triple mirrors

Reflect the light into itself over an angular range of about 12°, enabling extremely simple adjustment of the setup.

### (2.1) Large triple mirror in holder; on stand rod (9 mm x 10 mm dia.)

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for experiments on determining the velocity of light:

- 1 Lens,  $f = 200$  mm, dia. 120 mm (460 10)
- 1 Optical bench with standard profile, 0.5 m (from 460 34)
- 2 Optics riders, 60 mm high (460 351)
- 1 Ruler, e.g. Wooden ruler, 1 m (311 03) or Steel tape measure, 2 m (311 77)
- 1 Stand base, V-shape, large (300 01)
- 1 Stand rod, 100 cm (300 44)
- 1 Leybold-multiclamp (301 01)

or for measuring transit times on coaxial cables (see Fig. 4.2)

- 1 Coaxial cable with BNC-connectors, 50  $\Omega$ . 10 m long
- 1 T-connector, BNC (501 091)
- 1 Stand base, V-shape, small (300 02)

additionally required for reflection at a terminated end:

- 1 Straight, BNC (501 10)
- 1 Adapter BNC/4 mm, 2-pole (575 35)
- 1 Bridging plug (from 501 48)
- 1 STE-regulation resistor, 1 k $\Omega$  (577 86)



Fig. 5.2

Fig. 5.1

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Experiment setup

- Determining the velocity of light (Fig. 5.1); connection between 10-MHz output and second oscilloscope channel only required for measurements according to section 4.3.3
- Investigating transit times in coaxial cables (Fig. 5.2); for reflections at a terminated end, short the inner and outer conductors using a bridging plug (from 501 48) as shown in the detail drawing or connect them using a 1-k $\Omega$ -regulation resistor (577 86)

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## Notes on the listed equipment

### - Oscilloscope

The operating information given in this Instruction Sheet refers to the two-channel oscilloscope 303 (575 211) listed above. When using an oscilloscope with different technical specifications, you will need to adapt the settings and measuring instructions accordingly. The shape of the 10-MHz signals depends on the bandwidth of the oscilloscope (see section 3).

### - Connecting cables

The third BNC-cable is required for using the signals of the 10-MHz output on the second channel of the oscilloscope as the time base (see section 4.3.3). For distortion-free signal display, 50- $\Omega$  coaxial cables are required (use only the specified cables or those with comparable properties) Set up the experiment as shown in Fig. 4.1 or Fig. 4.2.

### Basic settings on the two-channel oscilloscope (575 211)

#### Mode:

- single-channel operation  
or, when using the 10-MHz output as a time base (see section 4.3.3)
- two-channel mode, alternating

channel I: DC, 5 mV/cm to 100 mV/cm, zero line at bottom edge of grid

channel II (two-channel mode only, for measuring method according to section 4.3.3): DC, 0.1 V/cm, zero line in middle of screen

Triggering: external, RF, + (rising edge),

Trigger level: automatic

Time-base sweep: 0.2 s/cm, cal, X-expansion x 1

Maximum intensity

## 4.2 Adjustment for measuring the velocity of light (Fig. 5.1)

### Coarse adjustment

- Set up the configuration so that the window (1.6) is about 20 cm in front of the lens (make sure the midpoints of the window and the lens are at the same height).
- Set up the large triple mirror at the appropriate distance  $s/2$  for the experiment
- $s/2$  for measuring procedure 4.3.1: about 10 m to 20 m  
4.3.2: about 5 m  
4.3.3: about 15 m
- Set up the mirror so that its center is at the level of the optical axis and its surface is aligned perpendicular to this axis.

### Fine adjustment, locating the signal

- Switch on the operating unit (plug in the supply unit) and the oscilloscope.
- To check the beam path: observe the triple mirror from a vantage point just above the housing and lens. If the mirror is not illuminated or is only red at the edges, vary the position of the optical bench by turning it slightly or adjusting it vertically with the adjusting screws until the center of the mirror is illuminated.
- Set the signal which then appears on the oscilloscope to the maximum possible amplitude by means of slight changes in the optical adjustment, particularly by carefully varying the position of the lens.
- Now switch the expansion of the x-axis to x 10 and bring the measuring signal into the desired position using the X-position knob (see section 4.3).

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Absolute measurement of the transit time  $t$  of light for a path  $s$  from the time interval between measuring signal  $U_1$  and reference signal  $U_0$

### 4.3 Measuring the velocity of light

#### 4.3.1 Measuring the transit time $t$ for a given path $s$ (absolute measurement with reference signal; see Fig. 6))

Position the measuring signal  $U_1$  generated according to section 4.2 in the right-hand area of the screen using the X-position knob.

To generate the reference signal  $U_0$ , place the small triple mirror on window (1.7); by moving the mirror on the opening, vary the effective reflecting surface so that  $U_0$  has exactly the same amplitude as  $U_1$ .

Use the X-position knob to shift the rising edge of reference signal  $U_0$  so that it intercepts the center line at its intersection with a grid line.

Read off transit time  $t$  at the point at which the measuring signal  $U_1$  intersects the center line.

*Note:*

The time difference between the reference pulse and the measuring pulse will only correspond to the difference between the two rising edges on the oscilloscope screen when both signals have the same amplitude and their separation is significantly greater than the pulse width.

#### 4.3.2 Measuring the change in transit time $\Delta t$ as a function of the change in path $s$ (relative measurement; see Fig. 7)

Position the measuring signal  $U_1$  generated according to section 4.2 using the X-position knob so that the maximum lies on a vertical grid line in the left-hand area of the screen (see Fig. 7.1).

After marking its position, move the triple mirror by  $\Delta s/2$ , so that the light path is increased by  $\Delta s$ . Determine the corresponding increase in the transit time  $\Delta t$  from the shift of the time signal (measurement in maximum, see Fig. 7.2)

Record multiple  $\Delta s$ - $\Delta t$ -value pairs and determine the velocity of light from the slope of the straight line  $\Delta t = f(\Delta s)$ .


Fig. 7.1/7.2

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Fig. 7.1/7.2

Relative measurement of the change in the transit time  $\Delta t$  of the light when changing the path by  $\Delta s$  from the change in the position of measuring signal  $U_1$

Fig. 7.1 Position of  $U_1$  at start of experiment

Fig. 7.2 Position of  $U_1$  after changing the path by  $\Delta s$


Fig. 8

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Absolute measurement of the velocity of light with the aid of an external time base  $t$  of 100 ns (period of the 10-MHz signal) from the light velocity measuring instrument.

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Display measuring signal  $U_1$  and reference signal  $U_0$  as described in section 4.3.1 and set them to exactly the same amplitude.

Additionally, connect the 10-MHz output and oscilloscope channel II according to section 4.1.

Using the phase adjuster (1.4), shift the 10-MHz signal so that the rising edge of the first pulse is on top of the rising edge of reference signal  $U_0$  as shown in Fig. 8.

Then, by moving the large triple mirror, vary the path  $s/2$  so that the rising edge of the measuring signal  $U_1$  precisely coincides with the second rising edge of the 10-MHz signal as shown in Fig. 8.

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#### 4.4 Measuring the transit time for signals on coaxial cables (see Fig. 5.2)

Set up the experiment as shown in Fig. 5.2. Via the T-connector (501 091), feed the signals generated with the small triple mirror to the oscilloscope as reference signal  $U_0$  and to the cable as the measuring pulse; this measuring pulse passes through the cable with length  $l$  twice, with reflection at the open or terminated end, and is displayed on the oscilloscope as measuring pulse  $U_1$  after transit time  $t$  has elapsed.

Determine the time  $t$  as described in section 4.3.1.