
Kontinuierliche Nebelkammer
Continuously Sensitive Cloud Chamber
Chambre à brouillard permanent

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1. Description

The continuously sensitive cloud chamber is used to observe the tracks of radioactive decay products with the aid of the condensation of cloud droplets on the ions formed on the tracks.

In contrast to the so-called expansion chamber (e. g. Cat. No. 559 57), this chamber is sensitive for several hours to the condensation on the tracks of the decay products. Therefore in continuous observation it is possible to render visible the cloud traces of even the weakest of radioactive samples.

In operation, the alcohol-moistened felt ring located under the cover constantly gives off vapour. By cooling the base, a temperature drop is maintained between the cover and the base which leads to the condensation of the falling vapour. In the stationary state, which develops after 20 to 30 minutes, there is a supersaturation of alcohol vapour a few centimetres above the base. Each ionising particle reaching this zone produces a cloud trace.

The chamber comprises the following parts:

Notes

1. The five-figure numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.
2. The data and illustrations are not binding in every detail for the design of the apparatuses. It is our sole aim always to adapt our manufacturing programme to the most recent knowledge gained in all scientific and technical fields.

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In addition we supply 2 preparations with the chamber:

- 1 needle with radium source for the simultaneous observation of alpha, beta and gamma rays.
1 needle coated with a very weak alpha-emitter.

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2. Operation

After being filled with dry ice, it takes the chamber about forty-five minutes before it becomes fully efficient. In order to save as much time as possible, the assembly etc. necessary for its operation is best carried out in the given sequence:

- Cover the base plate with a 1—2 mm layer of methanol.
- Moisten the felt ring of the cover plate with 100 ccs of a mixture prepared from 80 % methanol and 20 % water.
- Insert the silicon ring in the groove of the base plate, and the perbunan ring in the groove of the cover plate. Before clean the grooves and rings to get rid of any dust.
- Locate the glass cylinder. The light admission window of the glass cylinder should be placed between the locking screws of the light admission shaft. Close the lateral borings in the glass cylinder with rubber stoppers.
- Place on the cover plate and secure with three spring clips. On one side the space between two legs must be left free to locate the dry ice tray.
- The felt of the light admission shaft and the admission window to the glass cylinder should be

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well moistened with alcohol; place the light admission shaft in position so that it fits tightly.

- g) Take a block of dry ice approximately 180 mm in diameter and approx. 40 mm thick, and place it in the tray; press it down, and slide the tray under the chamber until firmly in position. If the ice is thicker, raise the chamber slightly. It is advisable to moisten the upper surface of the dry ice beforehand with alcohol.

h) Assembly of the illuminating system:

A source of light should be used consisting of a lamp housing (450 60a) with lamp (450 51) and a condenser made up of the lenses 75/150 (460 08) and 120/200 (460 10). The diaphragm should be placed on the lens (460 10). The light should come out of the condenser horizontally and weakly divergent, and pass through the light shaft into the chamber thus illuminating the section just above the bottom. These parts are best set up on a small optical bench (460 00) secured by means of a table clamp (301 06).

- i) To clean the chamber of ions, a D.C. of approximately 200 volts is required between the cover and the base. It is best obtained from the power unit (522 35) via a safety resistance (e. g. 536 25). To connect up the apparatus, there is a clamp on the side of the base and a boring on the side of the cover.

From time to time the stoppers in the lateral borings should be loosened to compensate for the low pressure in the chamber resulting from cooling. The intensive condensation shows that air has been drawn in when loosening the stoppers, and this is a reliable indication that the chamber is tight and that it is functioning properly. Should the chamber be leaky and not show a lower pressure, the rubber or silicon gaskets must be checked. Mostly it is found that such a condition is caused by a particle of dirt on one of the gaskets. In cases where difficulties are encountered in overcoming leakiness, it is advisable to rub the rubber gaskets with glycerine.

About half an hour after the cloud chamber has been switched on it commences to become sensitive. In a darkened room, with the illuminating system switched on, single traces of high-energy electrons appear which are caused by the cosmic radiation. Subsequently the electric field should

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Observations:

Once the chamber has become fully sensitive, paths may be constantly observed. The best direction of observation is against the light incidence through the wall of the glass cylinder. A plate with a weak radioactive source (e. g. 546 24) can be placed close to, and its secondary electrons observed. The needle with the radium source supplied with the apparatus should be held for a short period in front of a lateral boring of the glass cylinder opened for this purpose. After each of these experiments the electric field should be switched on, and the chamber cleaned of all ions occurring inside it. The weak alpha emitter supplied with the apparatus should be inserted in the chamber through one of the two lateral borings and moved constantly within the sensitive layer so that it is constantly in saturated vapour. The emission of the alpha-particles can be followed most clearly.

When inserting the preparation the chamber should not be open longer than absolutely necessary.

A little thorium emanation from the thorium source of the ionisation chamber (546 25) can be blown carefully into the chamber. In this experiment the electric cleansing field should be left connected to the chamber until the emanation has decreased sufficiently. The use of too much thorium emanation or of preparations which are too strong leads to a weakening of alcohol vapour in the sensitive layer so that the chamber is no longer able to indicate.

Further, it is possible to observe the secondary electrons emanating from X-rays. For this test, an X-ray tube (555 71) should be set up in front of the chamber and screened by means of a lead diaphragm with boring. When switched on, each time for a short period, the chamber is filled with the rippled paths of the secondary low-energy electrons.

All these, or similar experiments cannot, however, be carried out until the radiation background clearly visible in the chamber indicates the formation of the sensitive zone.

Dependent upon the outside conditions the dry ice tray must be refilled after 90 to 120 minutes. The block should not be allowed to become thinner than about 1 cm other-

wise the chamber will be too warm and will have to be cooled down to the normal pressure.

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wise the temperature of the base increases which means that the chamber would first have to be cooled again.

The chamber should be emptied on completion of the experiments, and left open until it is completely dry.

Special points:

Care should be taken not to contaminate the chamber even with weak activities since they can disturb or even completely ruin the proper functioning of the chamber. Above all radium emanation or similar long-lived substances should not be placed inside the chamber. When inserting the sources it is the edges of the openings in particular which are endangered. As a protective measure they are coated with nitro lacquer which, if necessary, can be scraped off, washed out with a cloth moistened with acetone and replaced by a new coating. This method should be applied in cases where an exceptionally strong activity, which is most discernible with alpha-radiation, occurs close to an opening. The outer finish of the glass cylinder should, however, be protected against acetone.

With respect to other countries, our agencies will be only too pleased to provide information concerning the possibilities of obtaining the dry ice required for the operation of the chamber (disk of 180 mm diameter).

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ASCA APARELHOS CIENTÍFICOS S.A.