5/1988 -Sf-

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The ionization chamber and telescopic cylinder are used to demonstrate and to examine radioactive materials. They can be used in particular to determine the ionizing effect of radioactive materials, the half-life of thorium emanation and the range of α particles.

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Literature:
Physics Experiments, Volume 3, "Atomic and Nuclear Physics" (599 942)
New Physics Leaflets for Colleges and Universities, Volume 1 "Nuclear Physics" (599 952)
Book "I-Measuring Amplifier D" (532 032)

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1 Safety Instructions

- Use a voltage which is not hazardous on contact (see para. 3)

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- ① Bottom of the ionization chamber with 4 mm socket for connecting the positive pole of the voltage supply and a circular insulator (inner diameter 13 mm), through which an electrode connected to the current sensitive measuring set up is led.
- ② Chamber lid with tube valves for injecting thoron (out of the thorium source 546 36) or gases upon which the ionization effects of radiation are to be examined.

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- ③ Metal pipe with sliding lid, together with bottom of chamber ① to make an ionization chamber of variable height

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Height: adjustable from 8.5 - 11 cm
Diameter: 7.5 cm

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3 Use

Important: a very sensitive ammeter is required due to the very low ionization current of 10^{-10} A to 10^{-11} A. Electrostatic induction may interfere noticeably.

Remedy:

Ensure proper earthing. It is helpful to connect the demonstration apparatus to the earthed socket using a second connecting rod and experimental cable.

Additionally required:

Contact safe voltage supply of approx. 300 V or 1 kV (adjustable for recording the characteristic curve); e.g.

300V power supply, 0.2 mA.....522 27
High voltage power supply, 10 kV.....522 37

Current sensitive measuring arrangement alternative

Wulf electroscop (Fig. 4).....546 00
or (Fig. 2)
I measuring amplifier D.....532 00
or (Fig. 3)
Electrometer amplifier.....532 14
Resistor, 10 G Ω577 03
Connecting rod.....532 16
Coupling plug.....340 89

Indicator or recorder for the amplifier, e.g.:

Moving coil measuring instrument D.....531 781
TY recorder.....575 70

Radioactive materials, e.g.:

Protected radium source with shield.....546 45
Am-241 preparation.....559 82
Thorium source.....546 36

Stands (see Fig.2/4 for examples)

Connect the ionization chamber to the voltage supply and the current sensitive measuring set up as shown in Fig.2/4.

When assembling the chamber on the electrometer amplifier, (see Fig.3), select the height of the bottom of the chamber on the connecting rod such that the upper surface of the chamber lid is at least 2 cm above the connecting rod or the attached preparation.

Fasten the radium preparation in the chamber (on a 4 mm pin) as shown in Figs. 2-4, using a coupling (e.g. from 501 641) if necessary.

Connect tube of the thorium source (546 36) with the tube valve as shown in Fig.3. Bring thoron into the chamber immediately before taking measurements by pressing the plastic bottle firmly several times.

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Recording the characteristic curve for the ionization chamber.

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Thoron half-life
Measurement of the ionization current against time using a moving coil instrument, measuring range 3 V DC (e.g. 531 781) and stop clock (e.g. 313 05)
or
display using a TY recorder (e.g. 575 70) connected to the output of an electrometer amplifier (532 14)

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